

# RTCB Antibody, HRP conjugated



PACO48195

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Catalytic subunit of the tRNA-splicing ligase complex that acts by directly joining spliced tRNA halves to mature-sized tRNAs by incorporating the precursor-derived splice junction phosphate into the mature tRNA as a canonical 3',5'-phosphodiester. May act as an RNA ligase with broad substrate specificity, and may function toward other RNAs.

**Gene ID:**

RTCB

**Uniprot**

Q9Y3I0

**Synonyms:**

tRNA-splicing ligase RtcB homolog (EC 6.5.1.3), RTCB, C22orf28

**Immunogen:**

Recombinant Human tRNA-splicing ligase RtcB homolog protein (52-229AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A