

LIPG Antibody, FITC conjugated



PACO48200

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Has phospholipase and triglyceride lipase activities. Hydrolyzes high density lipoproteins (HDL) more efficiently than other lipoproteins. Binds heparin.

Gene ID:

LIPG

Uniprot

Q9Y5X9

Synonyms:

Endothelial lipase (EC 3.1.1.3) (Endothelial cell-derived lipase) (EDL) (EL), LIPG

Immunogen:

Recombinant Human Endothelial lipase protein (360-500AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A