

GAPDH Antibody, HRP conjugated



PACO48203

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

GAPDH

Uniprot

P04406

Synonyms:

Glyceraldehyde-3-phosphate dehydrogenase (GAPDH) (EC 1.2.1.12) (Peptidyl-cysteine S-nitrosylase GAPDH) (EC 2.6.99), GAPDH, GAPD

Immunogen:

Recombinant Human Glyceraldehyde-3-phosphate dehydrogenase protein (2-335AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A