

ADORA2B Antibody, HRP conjugated



PACO48215

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Receptor for adenosine. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.

Gene ID:

ADORA2B

Uniprot

P29275

Synonyms:

Adenosine receptor A2b, ADORA2B

Immunogen:

Recombinant Human Adenosine receptor A2b protein (292-332AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A