

AGTR1 Antibody, HRP conjugated



PACO48223

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Receptor for angiotensin II. Mediates its action by association with G proteins that activate a phosphatidylinositol-calcium second messenger system.

Gene ID:

AGTR1

Uniprot

P30556

Synonyms:

Type-1 angiotensin II receptor (AT1AR) (AT1BR) (Angiotensin II type-1 receptor) (AT1), AGTR1, AGTR1A AGTR1B AT2R1 AT2R1B

Immunogen:

Recombinant Human Type-1 angiotensin II receptor protein (306-356AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A