

LIPE Antibody, HRP conjugated



PACO48383

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

In adipose tissue and heart, it primarily hydrolyzes stored triglycerides to free fatty acid, while in steroidogenic tissues, it principally converts cholesteryl esters to free cholesterol for steroid hormone production.

Gene ID:

LIPE

Uniprot

Q05469

Synonyms:

Hormone-sensitive lipase (HSL) (EC 3.1.1.79), LIPE

Immunogen:

Recombinant Human Hormone-sensitive lipase protein (1-251AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A