

# LIPE Antibody, FITC conjugated



PACO48384

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

In adipose tissue and heart, it primarily hydrolyzes stored triglycerides to free fatty acid, while in steroidogenic tissues, it principally converts cholesteryl esters to free cholesterol for steroid hormone production.

**Gene ID:**

LIPE

**Uniprot**

Q05469

**Synonyms:**

Hormone-sensitive lipase (HSL) (EC 3.1.1.79), LIPE

**Immunogen:**

Recombinant Human Hormone-sensitive lipase protein (1-251AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A