

# RPS6KA4 Antibody, FITC conjugated



PACO48456

---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	
<b>Reactivity:</b>	<b>Gene ID:</b>
Human	RPS6KA4
<b>Source:</b>	<b>Uniprot</b>
Rabbit	O75676
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	Ribosomal protein S6 kinase alpha-4 (S6K-alpha-4) (EC 2.7.11.1) (90 kDa ribosomal protein S6 kinase 4) (Nuclear mitogen- and stress-activated protein kinase 2) (Ribosomal protein kinase B) (RSKB), RPS6KA4, MSK2
<b>Applications:</b>	<b>Immunogen:</b>
ELISA	Recombinant Human Ribosomal protein S6 kinase alpha-4 protein (237-393AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A