

SFRP4 Antibody, FITC conjugated



PACO48464

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Soluble frizzled-related proteins (sFRPS) function as modulators of Wnt signaling through direct interaction with Wnts. They have a role in regulating cell growth and differentiation in specific cell types. SFRP4 may act as a regulator of adult uterine morphology and function. Increases apoptosis during ovulation possibly through modulation of FZ1/FZ4/WNT4 signaling. Has phosphaturic effects by specifically inhibiting sodium-dependent phosphate uptake.

Gene ID:

SFRP4

Uniprot

Q6FHJ7

Synonyms:

Secreted frizzled-related protein 4 (sFRP-4) (Frizzled protein, human endometrium) (FrpHE), SFRP4, FRPHE

Immunogen:

Recombinant Human Secreted frizzled-related protein 4 protein (293-346AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A