

# SLC6A3 Antibody, Biotin conjugated



PACO48473

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Amine transporter. Terminates the action of dopamine by its high affinity sodium-dependent reuptake into presynaptic terminals.

**Gene ID:**

SLC6A3

**Uniprot**

Q01959

**Synonyms:**

Sodium-dependent dopamine transporter (DA transporter) (DAT) (Solute carrier family 6 member 3), SLC6A3, DAT1

**Immunogen:**

Recombinant Human Sodium-dependent dopamine transporter protein (161-237AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A