

BOLA3 Antibody



PACO48630

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

BOLA3

Uniprot

Q53S33

Synonyms:

BolA-like protein 3, BOLA3


Immunogen:

Recombinant Human BolA-like protein 3 protein (1-79AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO48630 at dilution of 1:100.