

# NOTUM Antibody, HRP conjugated



PACO48659

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Carboxylesterase that acts as a key negative regulator of the Wnt signaling pathway by specifically mediating depalmitoylation of WNT proteins. Serine palmitoylation of WNT proteins is required for efficient binding to frizzled receptors.

**Gene ID:**

NOTUM

**Uniprot**

Q6P988

**Synonyms:**

Palmitoleoyl-protein carboxylesterase NOTUM (EC 3.1.1.98) (hNOTUM), NOTUM

**Immunogen:**

Recombinant Human Palmitoleoyl-protein carboxylesterase NOTUM protein (100-235AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A