

TRAPPC9 Antibody, HRP conjugated



PACO48703

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Functions as an activator of NF-kappa-B through increased phosphorylation of the IKK complex. May function in neuronal cells differentiation. May play a role in vesicular transport from endoplasmic reticulum to Golgi.

Gene ID:

TRAPPC9

Uniprot

Q96Q05

Synonyms:

Trafficking protein particle complex subunit 9 (NIK- and IKBKB-binding protein) (Tularik gene 1 protein), TRAPPC9, KIAA1882 NIBP

Immunogen:

Recombinant Human Trafficking protein particle complex subunit 9 protein (780-1030AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A