

# ESF1 Antibody



PACO48734

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May constitute a novel regulatory system for basal transcription. Negatively regulates ABT1.

**Gene ID:**

ESF1

**Uniprot**

Q9H501

**Synonyms:**

ESF1 homolog (ABT1-associated protein), ESF1, ABTAP C20orf6

**Immunogen:**

Recombinant Human ESF1 homolog protein (22-234AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A