

PUF60 Antibody, FITC conjugated



PACO48768

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PUF60

Uniprot

Q9UHX1

Synonyms:

Poly(U)-binding-splicing factor PUF60 (60 kDa poly(U)-binding-splicing factor) (FUSE-binding protein-interacting repressor) (FBP-interacting repressor) (Ro-binding protein 1) (RoBP1) (Siah-binding protein 1) (Siah-BP1), PUF60, FIR ROBPI SIAHBP1

Immunogen:

Recombinant Human Poly(U)-binding-splicing factor PUF60 protein (302-413AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A