

ZSWIM6 Antibody



PACO48774

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Gene ID:

ZSWIM6

Uniprot

Q9HCJ5

Synonyms:

Zinc finger SWIM domain-containing protein 6, ZSWIM6, KIAA1577


Immunogen:

Recombinant Human Zinc finger SWIM domain-containing protein 6 protein (430-583AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HepG2 cells using PACO48774 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO48774 at dilution of 1:100.