

# FAF1 Antibody, FITC conjugated



PACO48804

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Potentiates but cannot initiate FAS-induced apoptosis.

**Gene ID:**

FAF1

**Uniprot**

Q9UNN5

**Synonyms:**

FAS-associated factor 1 (hFAF1) (UBX domain-containing protein 12) (UBX domain-containing protein 3A), FAF1, UBXD12, UBXN3A

**Immunogen:**

Recombinant Human FAS-associated factor 1 protein (490-629AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A