

FAF1 Antibody, Biotin conjugated



PACO48805

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Potentiates but cannot initiate FAS-induced apoptosis.

Gene ID:

FAF1

Uniprot

Q9UNN5

Synonyms:

FAS-associated factor 1 (hFAF1) (UBX domain-containing protein 12) (UBX domain-containing protein 3A), FAF1, UBXD12, UBXN3A

Immunogen:

Recombinant Human FAS-associated factor 1 protein (490-629AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A