

ANXA4 Antibody



PACO48842

Product Information

Size:

50ug

Reactivity:

Dog

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

Gene ID:

ANXA4

Uniprot

P50994

Synonyms:

Annexin A4 (36 kDa zymogen granule membrane-associated protein) (ZAP36) (Annexin IV) (Annexin-4) (Lipocortin IV), ANXA4, ANX4

Immunogen:

Recombinant Dog Annexin A4 protein (2-319AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A