

DAZL Antibody, HRP conjugated



PACO48879

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

RNA-binding protein, which is essential for gametogenesis in both males and females. Plays a central role during spermatogenesis. Acts by binding to the 3'-UTR of mRNA, specifically recognizing GUU triplets, and thereby regulating the translation of key transcripts.

Gene ID:

DAZL

Uniprot

Q92904

Synonyms:

Deleted in azoospermia-like (DAZ homolog) (DAZ-like autosomal) (Deleted in azoospermia-like 1) (SPGY-like-autosomal), DAZL, DAZH DAZL1 DAZLA SPGYLA

Immunogen:

Recombinant Human Deleted in azoospermia-like protein (130-250AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A