

# E2F3 Antibody, HRP conjugated



PACO48887

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

E2F3

**Uniprot**

Q96AR0

**Synonyms:**

E2F3 protein

**Immunogen:**

Recombinant Human E2F3 protein (1-108AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A