

E2F3 Antibody, FITC conjugated



PACO48888

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

E2F3

Uniprot

Q96AR0

Synonyms:

E2F3 protein

Immunogen:

Recombinant Human E2F3 protein (1-108AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A