

IFNW1 Antibody



PACO48966

Product Information

Size:

50ug

Reactivity:

Bovine

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

extracellular space, cytokine activity, type I interferon receptor binding, adaptive immune response, B cell differentiation, B cell proliferation, cytokine-mediated signaling pathway, humoral immune response, innate immune response, natural killer cell activation involved in immune response.

Gene ID:

IFNW1

Uniprot

P07352

Synonyms:

Interferon omega-1 (IFN-omega-c1) (Interferon alpha-II-1), IFNW1, IFNW

Immunogen:

Recombinant Bovine Interferon omega-1 protein (24-195AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A