

# SHC1 Antibody, HRP conjugated



PACO49047

---

## Product Information

|                               |   |
|-------------------------------|---|
| <b>Size:</b>                  | <b>Protein Background:</b>  |
| 50ug                          |   |
| <b>Reactivity:</b>            | <b>Gene ID:</b>   |
| Human                         | SHC1  |
| <b>Source:</b>                | <b>Uniprot</b>  |
| Rabbit                        | P29353  |
| <b>Isotype:</b>               | <b>Synonyms:</b>  |
| IgG                           | SHC-transforming protein 1 (SHC-transforming protein 3) (SHC-transforming protein A) (Src homology 2 domain-containing-transforming protein C1) (SH2 domain protein C1), SHC1, SHC SHCA |
| <b>Applications:</b>          | <b>Immunogen:</b>   |
| ELISA                         | Recombinant Human SHC-transforming protein 1 protein (207-370AA).   |
| <b>Recommended dilutions:</b> | <b>Storage:</b>   |
|                               | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  |

## Product Images

---

N/A

N/A