

# SOX9 Antibody, Biotin conjugated



PACO49053

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Plays an important role in the normal skeletal development. May regulate the expression of other genes involved in chondrogenesis by acting as a transcription factor for these genes.

**Gene ID:**

SOX9

**Uniprot**

P48436

**Synonyms:**

Transcription factor SOX-9, SOX9

**Immunogen:**

Recombinant Human Transcription factor SOX-9 protein (296-501AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A