

TNFAIP3 Antibody



PACO49066

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:2000-1:10000,
IF:1:50-1:200

Protein Background:**Gene ID:**

TNFAIP3

Uniprot

P21580

Synonyms:

Tumor necrosis factor alpha-induced protein 3 (TNF alpha-induced protein 3) (EC 2.3.2) (EC 3.4.19.12) (OTU domain-containing protein 7C) (Putative DNA-binding protein A20) (Zinc finger protein A20) [Cleaved into: A20p50; A20p37], TNFAIP3, OTUD7C


Immunogen:

Recombinant Human Tumor necrosis factor alpha-induced protein 3 protein (97-183AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot. Positive WB detected in: RAW264.7 whole cell lysate, HeLa whole cell lysate, A549 whole cell lysate. All lanes: TNFAIP3 antibody at 2.7µg/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 90 kDa. Observed band size: 90 kDa.



Immunofluorescent analysis of HeLa cells using PACO49066 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).