

TP53I11 Antibody



PACO49070

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

negative regulation of cell proliferation, response to stress.

Gene ID:

TP53I11

Uniprot

O14683

Synonyms:

Tumor protein p53-inducible protein 11 (p53-induced gene 11 protein), TP53I11, PIG11


Immunogen:

Recombinant Human Tumor protein p53-inducible protein 11 protein (1-41AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HeLa cells using PACO49070 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).