

TP53I11 Antibody, FITC conjugated



PACO49072

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

negative regulation of cell proliferation, response to stress.

Gene ID:

TP53I11

Uniprot

O14683

Synonyms:

Tumor protein p53-inducible protein 11 (p53-induced gene 11 protein), TP53I11, PIG11

Immunogen:

Recombinant Human Tumor protein p53-inducible protein 11 protein (1-41AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A