

TP53I11 Antibody, Biotin conjugated



PACO49073

Product Information

Size:	Protein Background:
50ug	negative regulation of cell proliferation, response to stress.
Reactivity:	Gene ID:
Human	TP53I11
Source:	Uniprot
Rabbit	O14683
Isotype:	Synonyms:
IgG	Tumor protein p53-inducible protein 11 (p53-induced gene 11 protein), TP53I11, PIG11
Applications:	Immunogen:
ELISA	Recombinant Human Tumor protein p53-inducible protein 11 protein (1-41AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A