

TPH1 Antibody, HRP conjugated



PACO49075

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytosol, indolalkylamine biosynthetic process.

Gene ID:

TPH1

Uniprot

P17752

Synonyms:

Tryptophan 5-hydroxylase 1 (EC 1.14.16.4) (Tryptophan 5-monoxygenase 1), TPH1, TPH TPRH TRPH

Immunogen:

Recombinant Human Tryptophan 5-hydroxylase 1 protein (257-405AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A