

CDX2 Antibody, Biotin conjugated



PACO49117

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in the transcriptional regulation of multiple genes expressed in the intestinal epithelium. Important in broad range of functions from early differentiation to maintenance of the intestinal epithelial lining of both the small and large intestine.

Gene ID:

CDX2

Uniprot

Q99626

Synonyms:

Homeobox protein CDX-2 (CDX-3) (Caudal-type homeobox protein 2), CDX2, CDX3

Immunogen:

Recombinant Human Homeobox protein CDX-2 protein (12-116AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A