

# AQP1 Antibody, FITC conjugated



PACO49136

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Forms a water-specific channel that provides the plasma membranes of red cells and kidney proximal tubules with high permeability to water, thereby permitting water to move in the direction of an osmotic gradient.

**Gene ID:**

AQP1

**Uniprot**

P29972

**Synonyms:**

Aquaporin-1 (AQP-1) (Aquaporin-CHIP) (Urine water channel) (Water channel protein for red blood cells and kidney proximal tubule), AQP1, CHIP28

**Immunogen:**

Recombinant Human Aquaporin-1 protein (2-269AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A