

LOLPIB Antibody



PACO49194

Product Information

Size:

50ug

Reactivity:

Lolium perenne

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

LOLPIB

Uniprot

Q40240

Synonyms:

Major pollen allergen Lol p 5a (Allergen Lol p Ib) (Allergen Lol p Va) (allergen Lol p 5a), LOLPIB

Immunogen:

Recombinant Lolium perenne Major pollen allergen Lol p 5a protein (26-307AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A