

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probable substrate-recognition component of a SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins. Overexpression is leading to degradation of CBX5 and CBX1.

Gene ID:

Fbxw10

Uniprot

Q5SUS0

Synonyms:

F-box/WD repeat-containing protein 10 (F-box and WD-40 domain-containing protein 10), Fbxw10

Immunogen:

Recombinant Mouse F-box/WD repeat-containing protein 10 protein (701-1030AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A