

TP53INP1 Antibody, Biotin conjugated



PACO49241

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

TP53INP1

Uniprot

Q96A56

Synonyms:

Tumor protein p53-inducible nuclear protein 1 (Stress-induced protein) (p53-dependent damage-inducible nuclear protein 1) (p53DINP1), TP53INP1, P53DINP1 SIP

Immunogen:

Recombinant Human Tumor protein p53-inducible nuclear protein 1 protein (138-205AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A