

# UPF1 Antibody, Biotin conjugated



PACO49249

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

UPF1

**Uniprot**

Q92900

**Synonyms:**

Regulator of nonsense transcripts 1 (EC 3.6.4) (ATP-dependent helicase RENT1) (Nonsense mRNA reducing factor 1) (NORF1) (Up-frameshift suppressor 1 homolog (hUpf1), UPF1, KIAA0221 RENT1

**Immunogen:**

Recombinant Human Regulator of nonsense transcripts 1 protein (956-1125AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A