

STAC Antibody, HRP conjugated



PACO49275

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probably involved in a neuron-specific signal transduction.

Gene ID:

STAC

Uniprot

Q99469

Synonyms:

SH3 and cysteine-rich domain-containing protein (Src homology 3 and cysteine-rich domain-containing protein), STAC

Immunogen:

Recombinant Human SH3 and cysteine-rich domain-containing protein (1-119AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A