

HRH3 Antibody, Biotin conjugated



PACO49365

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

The H3 subclass of histamine receptors could mediate the histamine signals in CNS and peripheral nervous system. Signals through the inhibition of adenylate cyclase and displays high constitutive activity (spontaneous activity in the absence of agonist). Agonist stimulation of isoform 3 neither modified adenylate cyclase activity nor induced intracellular calcium mobilization.

Gene ID:

HRH3

Uniprot

Q9Y5N1

Synonyms:

Histamine H3 receptor (H3R) (HH3R) (G-protein coupled receptor 97), HRH3, GPCR97

Immunogen:

Recombinant Human Histamine H3 receptor protein (236-332AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A