

ABCC4 Antibody



PACO49366

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

May be an organic anion pump relevant to cellular detoxification.

Gene ID:

ABCC4

Uniprot

O15439

Synonyms:

Multidrug resistance-associated protein 4 (ATP-binding cassette sub-family C member 4) (MRP/cMOAT-related ABC transporter) (Multi-specific organic anion transporter B) (MOAT-B), ABCC4, MRP4


Immunogen:

Recombinant Human Multidrug resistance-associated protein 4 protein (562-697AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunofluorescent analysis of HepG2 cells using PACO49366 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).