

GCM2 Antibody, Biotin conjugated



PACO49629

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probable transcriptional regulator.

Gene ID:

GCM2

Uniprot

O75603

Synonyms:

Chorion-specific transcription factor GCMb (hGCMb) (GCM motif protein 2) (Glial cells missing homolog 2), GCM2, GCMB

Immunogen:

Recombinant Human Chorion-specific transcription factor GCMb protein (156-451AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A