

# IRF6 Antibody, Biotin conjugated



PACO49701

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Probable DNA-binding transcriptional activator. Key determinant of the keratinocyte proliferation-differentiation switch involved in appropriate epidermal development. Plays a role in regulating mammary epithelial cell proliferation. May regulate WDR65 transcription.

**Gene ID:**

IRF6

**Uniprot**

O14896

**Synonyms:**

Interferon regulatory factor 6 (IRF-6), IRF6

**Immunogen:**

Recombinant Human Interferon regulatory factor 6 protein (43-285AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A