

PACO49810

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Ras proteins bind GDP/GTP and possess intrinsic GTPase activity.

Gene ID:

NRAS

Uniprot

P01111

Synonyms:

GTPase NRas (Transforming protein N-Ras), NRAS, HRAS1


Immunogen:

Recombinant Human GTPase NRas protein (1-186AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO49810 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO49810 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).