

ONECUT2 Antibody, FITC conjugated



PACO49836

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transcriptional activator. Activates the transcription of a number of liver genes such as HNF3B.

Gene ID:

ONECUT2

Uniprot

O95948

Synonyms:

One cut domain family member 2 (Hepatocyte nuclear factor 6-beta) (HNF-6-beta) (One cut homeobox 2) (Transcription factor ONECUT-2) (OC-2), ONECUT2, HNF6B

Immunogen:

Recombinant Human One cut domain family member 2 protein (214-331AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A