

P2RX1 Antibody, FITC conjugated



PACO49848

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Ligand-gated ion channel with relatively high calcium permeability. Binding to ATP mediates synaptic transmission between neurons and from neurons to smooth muscle. Seems to be linked to apoptosis, by increasing the intracellular concentration of calcium in the presence of ATP, leading to programmed cell death.

Gene ID:

P2RX1

Uniprot

P51575

Synonyms:

P2X purinoceptor 1 (P2X1) (ATP receptor) (Purinergic receptor), P2RX1, P2X1

Immunogen:

Recombinant Human P2X purinoceptor 1 protein (54-175AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A