

# PLA2G10 Antibody, Biotin conjugated



PACO49877

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

PA2 catalyzes the calcium-dependent hydrolysis of the 2-acyl groups in 3-sn-phosphoglycerides. Has a powerful potency for releasing arachidonic acid, from cell membrane phospholipids. Prefers phosphatidylethanolamine and phosphatidylcholine liposomes to those of phosphatidylserine.

**Gene ID:**

PLA2G10

**Uniprot**

O15496

**Synonyms:**

Group 10 secretory phospholipase A2 (EC 3.1.1.4) (Group X secretory phospholipase A2) (GX sPLA2) (sPLA2-X) (Phosphatidylcholine 2-acylhydrolase 10), PLA2G10

**Immunogen:**

Recombinant Human Group 10 secretory phospholipase A2 protein (43-165AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A