

# PNLIPRP1 Antibody



PACO49886

---

## Product Information

**Size:**

50ug

**Reactivity:**

Dog

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May function as inhibitor of dietary triglyceride digestion. Lacks detectable lipase activity towards triglycerides, diglycerides, phosphatidylcholine, galactolipids or cholesterol esters (in vitro).

**Gene ID:**

PNLIPRP1

**Uniprot**

P06857

**Synonyms:**

Inactive pancreatic lipase-related protein 1 (PL-RP1), PNLIPRP1, PLRP1

**Immunogen:**

Recombinant Dog Inactive pancreatic lipase-related protein 1 protein (18-467AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A