

PRKDC Antibody, FITC conjugated



PACO49900

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PRKDC

Uniprot

P78527

Synonyms:

DNA-dependent protein kinase catalytic subunit (DNA-PK catalytic subunit) (DNA-PKcs) (EC 2.7.11.1) (DNPK1) (p460), PRKDC, HYRC HYRC1

Immunogen:

Recombinant Human DNA-dependent protein kinase catalytic subunit protein (3747-4015AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A