

PRL Antibody




PACO49902

Product Information

Size:	Protein Background:
50ug	Prolactin acts primarily on the mammary gland by promoting lactation.
Reactivity:	Gene ID:
Horse	PRL
Source:	Uniprot
Rabbit	P12420
Isotype:	Synonyms:
IgG	Prolactin (PRL), PRL
Applications:	Immunogen:
ELISA, WB	Recombinant Horse Prolactin protein (31-229AA).
Recommended dilutions:	Storage:
ELISA:1:2000-1:10000, WB:1:2000-1:10000	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: PRL antibody at 2.8 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. predicted band size: 28 kDa. observed band size: 28 kDa.