

PSMD9 Antibody, FITC conjugated



PACO49932

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Acts as a chaperone during the assembly of the 26S proteasome, specifically of the base subcomplex of the PA700/19S regulatory complex (RC). During the base subcomplex assembly is part of an intermediate PSMD9: PSMC6: PSMC3 module, also known as modulator trimer complex; PSMD9 is released during the further base assembly process.

Gene ID:

PSMD9

Uniprot

O00233

Synonyms:

26S proteasome non-ATPase regulatory subunit 9 (26S proteasome regulatory subunit p27), PSMD9

Immunogen:

Recombinant Human 26S proteasome non-ATPase regulatory subunit 9 protein (103-150AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A