

TERT Antibody, FITC conjugated



PACO50052

Product Information

Size:

50ug

Reactivity:

Dog

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

TERT

Uniprot

Q6A548

Synonyms:

Telomerase reverse transcriptase (EC 2.7.7.49) (Telomerase catalytic subunit), TERT

Immunogen:

Recombinant Dog Telomerase reverse transcriptase protein (1-308AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A